

## ELECTRODE PREPARATION GELS & PASTE



### ■ Pre-Tac

Skin preparation using this product increases conductivity and extends electrode wearing time. Helps prevent electrode movement under conditions of perspiration and humidity. Sold individually.

#### Product Description

#### Part Number

Pre-Tac 15ml bottle ..... 00103440



### ■ Prep-N-Stay

Skin preparation for the diaphoretic patient. Improves conductivity and prevents electrodes from falling off. Sold individually.

#### Product Description

#### Part Number

Prep-N-Stay 15ml bottle ..... 00103441

Prep-N-Stay 50ml bottle ..... 00103442



### ■ NuPrep

EEG and ECG skin Prep. Sold 3 per pack.

#### Product Description

#### Part Number

NuPrep 4 oz. Tube ..... 92100025



### ■ Lemon Prep

Lemon Prep is the only skin prep that combines the cleaning power of citrus and abrasion of pumice to remove oils and dead skin. This helps reduce impedances and improves electrode attachment and tracing. Sold individually.

#### Product Description

#### Part Number

Lemon Prep 4 oz. tube ..... 92000026



### ■ TEN20 Paste

EEG conductive paste that has adhesive qualities to hold electrodes in place. Sold 3 per pack.

#### Product Description

#### Part Number

TEN20 Paste 8 oz. Jar ..... 92100032

TEN20 Paste 4 oz. Jar ..... 92100030

TEN20 Paste 4 oz. Tube ..... 92100031



### ■ EC2 Cream

Conductive EEG electrode cream. Sold individually.

#### Product Description

#### Part Number

EC2 electrode cream 3.5 oz. Tube ... 92000000



### ■ Elefix Paste

Electrode adhesive, used with all surface applied electrodes. Sold individually.

Product Description	Part Number
Elefix 400g Tub .....	92000029
Elefix 180g Tube .....	92000028



### ■ TENSIVE® Conductive Adhesive Gel

Hypoallergenic, odorless conductive EEG adhesive gel.  
Sold individually.

Product Description	Part Number
TENSIVE 50g Tube .....	92000036



### ■ Signagel®

Highly conductive multi-purpose electrolyte gel. See thru tube lets you know exactly how much Signagel is left.  
Sold individually.

Product Description	Part Number
Signagel 8.5 oz. Tube .....	92000038



### ■ Quik-Gel Conductive Gel

Highly conductive electrolyte gel specially formulated for use with Compumedics Cap-based electrode placement systems.

Product Description	Part Number
Quik-Gel 16 oz. Jar .....	92000015
Quik-Gel 32 oz. Jar .....	92000016
Quik-Gel 128 oz. Jar .....	92000017



### ■ ECI Gel

Conductive electrolyte gel.

Product Description	Part Number
ECI Gel 16 oz. Jar .....	92000021
ECI Gel 32 oz. Jar .....	92000019
ECI Gel 128 oz. Jar .....	92000020